





Growable seat shell base frame



Tiltable wheelchair for individual positioning



# Loop SORG

#### The seat shell base

German aid index N° 26.99.01.1078



back height
43 or 50 cm



seat depth from 32 cm (+6 cm)



weight capacity
120 kg



wheel sizes 12", 20" - 24"



frame width

Size 1 30, 34, 38, 42cm

Size 2 34, 38, 42, 46cm

Size 3 38, 42, 46, 50cm




tilting angle -5° till +35° or +2,5° till + 40°



#### Overwiev

- Seat shell base frame that grows with the patient
- Tiltable from -5° to +35° or from +2,5° to +40°
- Parameters such as seat height, seat position/centre of gravity, distance of foot plate/s etc. can be adjusted precisely to the user
- alternatively available with 12" wheels
- suitable as a seat for transport in a vehicle according to ISO 7176-19





## Loop SORG RS

### The tiltable wheelchair

for your individual body guidance



seat width **30-48 cm** (+ 4 cm)



back height **30-60 cm** (+5 cm)



seat depth 30-48 cm (max. +6 cm)



weight capacity
120 kg



wheel sizes 20"-24"



tilting angle
-3° to max. +34°



#### Overview

- Tiltable active wheelchairs
- The optionally available ERGO seat and back system offers an individual and best possible body guidance
- The 3-part ERGO back is growing in the width and offers space for a seat designed individually by a specialist shop
- suitable as a seat for transport in a vehicle according to ISO 7176-19



### One name - two models

### Family overview:

|  | Loop <sup>sorg</sup> | Loop <sup>sorg</sup> RS |
|--|----------------------|-------------------------|
| seat depth   | from 32 cm           | 30 - 48 cm              |
| back height  | 43 or 50 cm          | 30 - 60 cm              |
| max. weight capacity   | 120 kg               | 120 kg                  |
| wheel sizes  | 12", 20" - 24"       | 20" - 24"               |
| caster sizes   | 5", 5,5", 6", 7"     | 5", 5,5", 6", 7"        |
| tilting angle  | ✓                    | <b>⊘</b>                |
| growable in seat width, seat height<br>and back height                 | <b>⊘</b>             | <b>⊘</b>                |
| suitable as a seat for transport in a vehicle according to ISO 7176-19 | <b>⊘</b>             | ✓                       |
| specially adapted seat shell from the specialized supplier necessary   | <b>⊘</b>             |                         |

## Equipment examples

All equipment options are on the order form of the corresponding wheelchair or sheat shell base



12" wheel



detachable leg supports



seat shell adapter with wedge (only for Loop<sup>SORG</sup>)



drum brake



back can be folded down to the front

Angl

much

mor3 ...



Outdoor Front End



ERGO-back



Push bail with synthetic hinge



ERGO-side guards with upholstery

