





Active, tiltable and growable





## Tilty Vario

## Active, tiltable, growable

German aid index N° 18.99.02.1013



seat width **28-44 cm** (+2 cm)



back height **30-50 cm** (+5 cm)



seat depth **30-44 cm** (+4 cm)



weight capacity

100 kg



wheel sizes 22" and 24"



tilting angle from 0° to +25°









## Overview

- · Tiltable wheelchair for adolescents and adults
- Back adjustable in angle (80° 120°) and completely stowable
- Highly stable rigid frame with 10 cm abduction; also suitable for Pörnbacher therapy methods
- Smooth handling with minimized effort
- Suitable as a seat for transport in a vehicle according to ISO 7176-19

## Equipment examples

All equipment options are on the order form of the Tilty Vario



holder to attach abduction wedge



truss pads



arm pads adjustable in height and angle



undivided foot rest



head rests with chestshoulder-belt

And much more ...